

Product datasheet for **GA110287**

Human TRPM7 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TRPM7 CRISPRa kit - CRISPR gene activation of human transient receptor potential cation channel subfamily M member 7
Symbol:	TRPM7
UniProt ID:	<u>Q96QT4</u>
Synonyms:	ALSPDC; CHAK; CHAK1; LTrpC-7; LTRPC7; TRP-PLIK
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)