


Product datasheet for **GA110290**

Human GIN1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	GIN1 CRISPRa kit - CRISPR gene activation of human gypsy retrotransposon integrase 1
Symbol:	GIN1
Components:	GA110290G1, GIN1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110290G2, GIN1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110290G3, GIN1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NXP7</u>
Synonyms:	GIN-1; TGIN1; ZH2C2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)