


Product datasheet for **GA110311**

Human PAQR5 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PAQR5 CRISPRa kit - CRISPR gene activation of human progesterin and adipoQ receptor family member 5
Symbol:	PAQR5
Components:	GA110311G1, PAQR5 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110311G2, PAQR5 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110311G3, PAQR5 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NXX6</u>
Synonyms:	MPRG
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)