

Product datasheet for **GA110312**

Human WDR55 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	WDR55 CRISPRa kit - CRISPR gene activation of human WD repeat domain 55
Symbol:	WDR55
Components:	GA110312G1, WDR55 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110312G2, WDR55 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110312G3, WDR55 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9H6Y2</u>
Synonyms:	FLJ20195; FLJ21702
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)