

Product datasheet for **GA110315**

Human GON4L activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	GON4L CRISPRa kit - CRISPR gene activation of human gon-4 like
Symbol:	GON4L
Components:	GA110315G1, GON4L gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110315G2, GON4L gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110315G3, GON4L gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q3T8J9
Synonyms:	GON-4; GON4; YARP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)