

Product datasheet for **GA110325**

Human TMEM214 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM214 CRISPRa kit - CRISPR gene activation of human transmembrane protein 214
Symbol:	TMEM214
Components:	GA110325G1, TMEM214 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110325G2, TMEM214 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110325G3, TMEM214 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6NUQ4</u>
Synonyms:	FLJ20254; FLJ39682
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)