

Product datasheet for **GA110329**

Human PIGG activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PIGG CRISPRa kit - CRISPR gene activation of human phosphatidylinositol glycan anchor biosynthesis class G
Symbol:	PIGG
UniProt ID:	Q5H8A4
Synonyms:	GPI7; LAS21; MRT53; PRO4405; RLGS1930
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)