

Product datasheet for **GA110338**

Human TEX10 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TEX10 CRISPRa kit - CRISPR gene activation of human testis expressed 10
Symbol:	TEX10
Components:	GA110338G1, TEX10 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110338G2, TEX10 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110338G3, TEX10 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NXF1</u>
Synonyms:	bA208F1.2; lpi1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)