

Product datasheet for **GA110340**

Human CWC25 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CWC25 CRISPRa kit - CRISPR gene activation of human CWC25 spliceosome associated protein homolog
Symbol:	CWC25
Components:	GA110340G1, CWC25 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110340G2, CWC25 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110340G3, CWC25 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NXE8</u>
Synonyms:	CCDC49
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

