

Product datasheet for **GA110341**

Human Retinol Saturase (RETSAT) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	RETSAT CRISPRa kit - CRISPR gene activation of human retinol saturase
Symbol:	RETSAT
UniProt ID:	<u>Q6NUM9</u>
Synonyms:	FLJ20296
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)