

Product datasheet for **GA110343**

Human plasticity related gene 3 (PLPPR1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLPPR1 CRISPRa kit - CRISPR gene activation of human phospholipid phosphatase related 1
Symbol:	PLPPR1
UniProt ID:	Q8TBJ4
Synonyms:	LPPR1; PRG-3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)