

Product datasheet for **GA110347**

Human INO80D activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	INO80D CRISPRa kit - CRISPR gene activation of human INO80 complex subunit D
Symbol:	INO80D
Components:	GA110347G1, INO80D gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110347G2, INO80D gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110347G3, INO80D gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q53TQ3</u>
Synonyms:	FLJ20309
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)