

Product datasheet for **GA110350**

Human RNF43 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	RNF43 CRISPRa kit - CRISPR gene activation of human ring finger protein 43
Symbol:	RNF43
UniProt ID:	Q68DV7
Synonyms:	RNF124; SSPCS; URCC
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)