


Product datasheet for **GA110370**

Human DUS2L (DUS2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DUS2 CRISPRa kit - CRISPR gene activation of human dihydrouridine synthase 2
Symbol:	DUS2
Components:	GA110370G1, DUS2L gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110370G2, DUS2L gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110370G3, DUS2L gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NX74</u>
Synonyms:	DUS2L; SMM1; URLC8
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)