


Product datasheet for **GA110378**

Human TMEM161A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM161A CRISPRa kit - CRISPR gene activation of human transmembrane protein 161A
Symbol:	TMEM161A
Components:	GA110378G1, TMEM161A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110378G2, TMEM161A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110378G3, TMEM161A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NX61</u>
Synonyms:	AROS-29
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)