

Product datasheet for **GA110384**

Human DUSP23 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DUSP23 CRISPRa kit - CRISPR gene activation of human dual specificity phosphatase 23
Symbol:	DUSP23
UniProt ID:	Q9BVJ7
Synonyms:	DUSP25; LDP-3; MOSP; VHZ
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)