

Product datasheet for **GA110403**

Human TXNL4B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TXNL4B CRISPRa kit - CRISPR gene activation of human thioredoxin like 4B
Symbol:	TXNL4B
Components:	GA110403G1, TXNL4B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110403G2, TXNL4B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110403G3, TXNL4B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NX01</u>
Synonyms:	Dim2; DLP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)