

Product datasheet for **GA110407**

Human SSH3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SSH3 CRISPRa kit - CRISPR gene activation of human slingshot protein phosphatase 3
Symbol:	SSH3
UniProt ID:	<u>Q8TE77</u>
Synonyms:	SSH3L
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)