


Product datasheet for **GA110408**

Human TIPIN activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TIPIN CRISPRa kit - CRISPR gene activation of human TIMELESS interacting protein
Symbol:	TIPIN
Components:	GA110408G1, TIPIN gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110408G2, TIPIN gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110408G3, TIPIN gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9BVW5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)