

Product datasheet for **GA110421**

Human SLC25A38 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SLC25A38 CRISPRa kit - CRISPR gene activation of human solute carrier family 25 member 38
Symbol:	SLC25A38
Components:	GA110421G1, SLC25A38 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110421G2, SLC25A38 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110421G3, SLC25A38 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q96DW6
Synonyms:	SIDBA2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)