

Product datasheet for **GA110426**

Human CLN6 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CLN6 CRISPRa kit - CRISPR gene activation of human CLN6 transmembrane ER protein
Symbol:	CLN6
Components:	GA110426G1, CLN6 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110426G2, CLN6 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110426G3, CLN6 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NWW5</u>
Synonyms:	CLN4A; HsT18960; nclf
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)