

## Product datasheet for **GA110463**

### Human HEATR3 activation kit by CRISPRa

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                               |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| <b>Description:</b>  | HEATR3 CRISPRa kit – CRISPR gene activation of human HEAT repeat containing 3                                                                                                                                                                                                                                 |
| <b>Symbol:</b>       | HEATR3                                                                                                                                                                                                                                                                                                        |
| <b>Components:</b>   | GA110463G1, HEATR3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA110463G2, HEATR3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA110463G3, HEATR3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| <b>UniProt ID:</b>   | <u><a href="#">Q7Z4Q2</a></u>                                                                                                                                                                                                                                                                                 |
| <b>Synonyms:</b>     | SYO1                                                                                                                                                                                                                                                                                                          |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:

