


Product datasheet for **GA110467**

Human SLC35A5 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SLC35A5 CRISPRa kit - CRISPR gene activation of human solute carrier family 35 member A5
Symbol:	SLC35A5
Components:	GA110467G1, SLC35A5 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110467G2, SLC35A5 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110467G3, SLC35A5 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9BS91
Synonyms:	DKFZp434E102; FLJ11130; FLJ20730; FLJ25973
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)