

## Product datasheet for **GA110472**

### Human PTC3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PTCD3 CRISPRa kit - CRISPR gene activation of human pentatricopeptide repeat domain 3
Symbol:	PTCD3
UniProt ID:	<a href="#">Q96EY7</a>
Synonyms:	MRP-S39
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)