

## Product datasheet for **GA110487**

### Human C9orf68 (SPATA6L) activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                  |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                     |
| Description:  | SPATA6L CRISPRa kit - CRISPR gene activation of human spermatogenesis associated 6 like                                                                                                                                                                                                                          |
| Symbol:       | SPATA6L                                                                                                                                                                                                                                                                                                          |
| Components:   | GA110487G1, C9orf68 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA110487G2, C9orf68 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA110487G3, C9orf68 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u>Q8N4H0</u>                                                                                                                                                                                                                                                                                                    |
| Synonyms:     | bA6J24.2; C9orf68                                                                                                                                                                                                                                                                                                |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                          |

#### Product images:



[View online »](#)