

Product datasheet for **GA110493**

Human C9orf40 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	C9orf40 CRISPRa kit - CRISPR gene activation of human chromosome 9 open reading frame 40
Symbol:	C9orf40
Components:	GA110493G1, C9orf40 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110493G2, C9orf40 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110493G3, C9orf40 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8IXQ3</u>
Synonyms:	FLJ10110; FLJ25795
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)