

Product datasheet for **GA110506**

Human SLC38A4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SLC38A4 CRISPRa kit - CRISPR gene activation of human solute carrier family 38 member 4
Symbol:	SLC38A4
Components:	GA110506G1, SLC38A4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110506G2, SLC38A4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110506G3, SLC38A4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q969I6
Synonyms:	ATA3; NAT3; PAAT; SNAT4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)