


Product datasheet for **GA110524**

Human **DYNC2I1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	WDR60 CRISPRa kit - CRISPR gene activation of human WD repeat domain 60
Symbol:	DYNC2I1
Components:	GA110524G1, DYNC2I1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110524G2, DYNC2I1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110524G3, DYNC2I1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8WVS4
Synonyms:	FAP163; SRPS6; SRTD8
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)