

Product datasheet for **GA110533**

Human PIWIL2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PIWIL2 CRISPRa kit - CRISPR gene activation of human piwi like RNA-mediated gene silencing 2
Symbol:	PIWIL2
UniProt ID:	Q8TC59
Synonyms:	CT80; HILI; mili; PIWIL1L
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)