

Product datasheet for **GA110542**

Human WDR79 (WRAP53) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	WRAP53 CRISPRa kit - CRISPR gene activation of human WD repeat containing antisense to TP53
Symbol:	WRAP53
UniProt ID:	Q9BUR4
Synonyms:	DKCB3; TCAB1; WDR79
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)