


Product datasheet for **GA110556**

Human TMEM38B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM38B CRISPRa kit - CRISPR gene activation of human transmembrane protein 38B
Symbol:	TMEM38B
Components:	GA110556G1, TMEM38B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110556G2, TMEM38B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110556G3, TMEM38B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NVV0</u>
Synonyms:	bA219P18.1; C9orf87; D4Ertd89e; OI14; TRIC-B; TRICB
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)