

Product datasheet for **GA110564**

Human TMEM33 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | TMEM33 CRISPRa kit - CRISPR gene activation of human transmembrane protein 33 |
| Symbol: | TMEM33 |
| UniProt ID: | P57088 |
| Synonyms: | 1600019D15Rik; SHINC-3; SHINC3 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)