

Product datasheet for **GA110586**

Human **SLC25A36** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SLC25A36 CRISPRa kit - CRISPR gene activation of human solute carrier family 25 member 36
Symbol:	SLC25A36
UniProt ID:	Q96CQ1
Synonyms:	PNC2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)