

Product datasheet for **GA110596**

Human P15RS (RPRD1A) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	RPRD1A CRISPRa kit - CRISPR gene activation of human regulation of nuclear pre-mRNA domain containing 1A
Symbol:	RPRD1A
Components:	GA110596G1, P15RS gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110596G2, P15RS gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110596G3, P15RS gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96P16</u>
Synonyms:	HsT3101; P15RS
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)