


## Product datasheet for **GA110599**

### Human MyoGEF (PLEKHG6) activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | PLEKHG6 CRISPRa kit - CRISPR gene activation of human pleckstrin homology and RhoGEF domain containing G6                                                                                                                                                                                                     |
| Symbol:       | PLEKHG6                                                                                                                                                                                                                                                                                                       |
| Components:   | GA110599G1, MyoGEF gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA110599G2, MyoGEF gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA110599G3, MyoGEF gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q3KR16</a>                                                                                                                                                                                                                                                                                        |
| Synonyms:     | MyoGEF                                                                                                                                                                                                                                                                                                        |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)