

Product datasheet for **GA110608**

Human ATAD3A activation kit by CRISPRa

Product data:

| | |
|---------------|---------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits |
| Description: | ATAD3A CRISPRa kit - CRISPR gene activation of human ATPase family AAA domain containing 3A |
| Symbol: | ATAD3A |
| UniProt ID: | Q9NVI7 |
| Synonyms: | FLJ10709; FLJ35514 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)