

Product datasheet for **GA110624**

Human Lano (LRRC1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LRRC1 CRISPRa kit - CRISPR gene activation of human leucine rich repeat containing 1
Symbol:	LRRC1
Components:	GA110624G1, Lano gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110624G2, Lano gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110624G3, Lano gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9BTT6</u>
Synonyms:	dj523E19.1; LANO
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)