


Product datasheet for **GA110637**

Human **SLC47A1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SLC47A1 CRISPRa kit - CRISPR gene activation of human solute carrier family 47 member 1
Symbol:	SLC47A1
Components:	GA110637G1, SLC47A1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110637G2, SLC47A1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110637G3, SLC47A1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96FL8</u>
Synonyms:	MATE1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)