

Product datasheet for **GA110647**

Human TMEM39A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM39A CRISPRa kit - CRISPR gene activation of human transmembrane protein 39A
Symbol:	TMEM39A
Components:	GA110647G1, TMEM39A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110647G2, TMEM39A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110647G3, TMEM39A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NV64</u>
Synonyms:	FLJ10902
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)