

Product datasheet for **GA110648**

Human **WDR41** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	WDR41 CRISPRa kit – CRISPR gene activation of human WD repeat domain 41
Symbol:	WDR41
Components:	GA110648G1, WDR41 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110648G2, WDR41 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110648G3, WDR41 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9HAD4</u>
Synonyms:	MSTP048
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

