

Product datasheet for **GA110679**

Human SAPS3 (PPP6R3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PPP6R3 CRISPRa kit - CRISPR gene activation of human protein phosphatase 6 regulatory subunit 3
Symbol:	PPP6R3
UniProt ID:	Q5H9R7
Synonyms:	C11orf23; PP6R3; SAP190; SAPL; SAPLa; SAPS3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)