

Product datasheet for **GA110687**

Human PI 4 Kinase type 2 beta (PI4K2B) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PI4K2B CRISPRa kit - CRISPR gene activation of human phosphatidylinositol 4-kinase type 2 beta
Symbol:	PI4K2B
UniProt ID:	<u>Q8TCG2</u>
Synonyms:	PI4KIIIB; PIK42B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)