

Product datasheet for **GA110690**

Human SPTLC3 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | SPTLC3 CRISPRa kit - CRISPR gene activation of human serine palmitoyltransferase long chain base subunit 3 |
| Symbol: | SPTLC3 |
| UniProt ID: | Q9NUV7 |
| Synonyms: | C20orf38; dj718P11; dj718P11.1; hLCB2b; LCB2B; LCB3; SPT3; SPTLC2L |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)