

Product datasheet for **GA110692**

Human ZNF444 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF444 CRISPRa kit - CRISPR gene activation of human zinc finger protein 444
Symbol:	ZNF444
Components:	GA110692G1, ZNF444 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110692G2, ZNF444 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110692G3, ZNF444 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8N0Y2</u>
Synonyms:	EZF-2; EZF2; ZSCAN17
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)