


Product datasheet for **GA110699**

Human TMA16 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMA16 CRISPRa kit - CRISPR gene activation of human translation machinery associated 16 homolog
Symbol:	TMA16
Components:	GA110699G1, TMA16 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110699G2, TMA16 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110699G3, TMA16 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96EY4</u>
Synonyms:	C4orf43
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)