

## Product datasheet for **GA110715**

### Human NIPSNAP3B activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                        |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                           |
| Description:  | NIPSNAP3B CRISPRa kit - CRISPR gene activation of human nipsnap homolog 3B                                                                                                                                                                                                                                             |
| Symbol:       | NIPSNAP3B                                                                                                                                                                                                                                                                                                              |
| Components:   | GA110715G1, NIPSNAP3B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA110715G2, NIPSNAP3B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA110715G3, NIPSNAP3B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q9BS92</a></u>                                                                                                                                                                                                                                                                                          |
| Synonyms:     | FP944; NIPSNAP3; SNAP1                                                                                                                                                                                                                                                                                                 |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                                |

#### Product images:



[View online »](#)