

Product datasheet for **GA110722**

Human SLC35C1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SLC35C1 CRISPRa kit - CRISPR gene activation of human solute carrier family 35 member C1
Symbol:	SLC35C1
Components:	GA110722G1, SLC35C1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110722G2, SLC35C1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110722G3, SLC35C1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q96A29
Synonyms:	CDG2C; FUCT1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)