

Product datasheet for **GA110736**

Human PI4K2A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PI4K2A CRISPRa kit - CRISPR gene activation of human phosphatidylinositol 4-kinase type 2 alpha
Symbol:	PI4K2A
UniProt ID:	<u>Q9BTU6</u>
Synonyms:	PI4KII; PIK42A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)