

Product datasheet for **GA110740**

Human TMEM176A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM176A CRISPRa kit - CRISPR gene activation of human transmembrane protein 176A
Symbol:	TMEM176A
Components:	GA110740G1, TMEM176A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA110740G2, TMEM176A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA110740G3, TMEM176A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q96HP8
Synonyms:	GS188; HCA112; MS4B1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)