

## Product datasheet for **GA110741**

### Human GPR48 (LGR4) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	LGR4 CRISPRa kit - CRISPR gene activation of human leucine rich repeat containing G protein-coupled receptor 4
Symbol:	LGR4
UniProt ID:	<a href="#">Q9BXB1</a>
Synonyms:	BNMD17; GPR48
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)